

VED-1
March 1, 1992

ACTION PLAN 1992

1. IMMUNOSCREENING OF EXPRESSION LIBRARY

2. Use amino acid sequence information to synthesize primers for use in PCR.

3. USE PCR to isolate desired sequences from the tobacco root RNA and DNA.

a. RNA PCR (Perkin-Elmer kit, Forward and Reverse primers)

b. DNA PCR (Perkin-Elmer kit, Forward and Reverse Primers)

c. ANCHORED PCR: AMPLIFICATION WITH SINGLE-SIDED SPECIFICITY

1. Take tobacco root poly A RNA. Synthesize complementary strand of DNA using oligo dT as primer. Destroy RNA by alkali hydrolysis. Amplify Using

2023141738

combinations of forward primers and oligo dT25 .CLONE the PCR PRODUCT.SEQUENCE the insert. Enter the sequence into computer.use DNA STAR(findpro programme) to find PMT1.PRO,PMT2.PRO, PMT3.PRO PMT4.PRO sequences.

2.ANCHORED PCR of cDNA with OLIGO(dG)tails.Synthesise cDNA using REVERSE PCR PRIMERS or oligo dT.Add oligo dG to the newly synthesised c DNA strand.Amplify using REVERSE PCR PRIMERS and OLIGO dC.Clone the PCR product AND analyse as in 1.

d.LIGATION-MEDIATED ANCHORED PCR.Digest the genomic DNA with a Mbo1 OR its isoschizomer Sau3A.

If possible, Fill one base to prevent self-Ligation.

Ligate an anchor-primer on the ends of the DNA fragments.

Amplify with a specific primer as well as the anchor.

The choice of which strand the primer is derived from determines which side of the fragment is amplified .

This can be done for several enzymes that leave cohesive ends with a four-base 5'overhang. This could allow isolation of the promoter for the gene of interest.

4. GENERATION OF SINGLE-STRANDED LABELED PRIMERS AND THEIR USE IN NORTHERN AND SOUTHERN HYBRIDIZATION AND SCREENING OF LIBRARIES.

5. PMT POLYCLONALS(GPC FRACTION)

IMMUNOSCREENING OF cDNA EXPRESSION LIBRARY.

1. EXPRESSION LIBRARY OF TOBACCO ROOT mRNA IS AT HAND.
2. A POLYCLONAL ANTIBODY AGAINST A SYNTHETIC PEPTIDE N-29 MER HAS BEEN RAISED
3. THIS POLYCLONAL ANTIBODY LIGHTS UP 60 KD REGION IN WESTERN BLOTS OF TOBACCO PROTEIN.
4. THE EXPRESSION LIBRARY(1) WILL BE SCREENED USING POLYCLONAL ANTIBODY(2).
5. IF A POSITIVE CLONE IS FOUND, THE PROTEIN OF THAT CLONE MAY BE PURIFIED USING ANTI betaGAL ANTIBODIES.

SEQUENCE INFORMATION

1. N-terminus of intact 60 kD protein, electroblot/sequenced. (code:PMT1.PRO)

PMT L S S N F L F G T A S S Y Y Q Y E G A F L S D G V G L S N --- COOH
Deg 6 6 6 2 2 6 2 4 4 4 6 6 2 2 2 2 2 4 4 2 6 6 2 4 4 4 6 6 2

Alignments:

PNMT(PET) S(16) S S Y Q
PEMT(PhosphatidyleMT) F Y G D F L S

2. 6.2 kD CNBr Fragment of 60 kD protein, electroblot/sequenced (F-SG-T MOTIF, region I), (a mat would precede all of the CNBr fragments). (code:PMT2.PRO)

PMT N E P X E V A I S G Y R D X T --- COOH
Deg 2 2 4 2 4 4 3 6 4 2 6 2 4

Alignments:

erm C/G
RNA MT I E E I G S G K G H F T
PNMT (65) E A T G E S G R T
PNMT (75) L I D I G S G P T Y N L

3. 3 kD CNBr Fragment of 60 kD after removal of N-terminus (DAL - Region II). (code:PMT3.PRO)

PMT A D I E H Y S K L I D A L X I K G I Q F --- COOH
G L L A P
Deg 2 6 2 6 3 2 4 6 3 2 4 3

Alignments:

PNMT (125) Y S K G
PNMT (174) P A D A L
AdoHcyHMT K D A I

4. 15 kD CNBr Fragment of 60 kD protein. (code:PMT4.PRO)

PMT F I T E N G F A G R S G R P (?) --- COOH
Deg 2 3 4 2 2 4 2 4 4 6 6 4 6

5. Possible sequence order of fragments obtained for 60 kD protein.

From this speculative "motif" examination, a deduced composite can be drawn for putative PMT as follows (cf. PNMT, PNAS 83,5455 (1986):

N-terminus (~3 kD) ----- 6.2 kD fragment ----- other 3 kD fragment ----- 15 kD fragment ----- COOH

Ref. Ingrosso et al., JBC 264,20131-20139 (1989); Baetge et al., PNAS 83,5454-5458. Also see. Lauster et al., FEBS Lett. 220, 167-176.

2023141742

PCR AND PRIMERS

by VED MALIK

1. Importance of first PCR cycle AND 3' end of primers
2. Important to know where you are
3. Examination of PMT. SEquences.
4. N-terminal codons included in primers
5. C-terminal codons should appear in amplified product.
6. C-terminal codons can constitute the primers, and N-terminal codons can appear in the PCR product.

2023141743

PMT1.PRO

VED-7

-- PAGE 1

PROBE#: 1

Melting Temp: 52

Number of possibilities: 128

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala 19

Oligonucleotide sequence: taytaycartaygargngc

PROBE#: 2

Melting Temp: 50

Number of possibilities: 128

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala 19

Oligonucleotide sequence: aytaycartaygargngc

PROBE#: 3

Melting Temp: 52

Number of possibilities: 256

PMT1.PRO 14 Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: taycartaygargngc

PROBE#: 4

Melting Temp: 50

Number of possibilities: 256

PMT1.PRO 14 Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: taycartaygargngc

PROBE#: 5

Melting Temp: 50

Number of possibilities: 256

PMT1.PRO 14 Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: aycartaygargngc

✓ PROBE#: 6

Melting Temp: 58

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: taytaycartaygargngc

PROBE#: 7

Melting Temp: 56

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: taytaycartaygargngc

PROBE#: 8

Melting Temp: 56

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: aytaycartaygargngc

PROBE#: 9

Melting Temp: 54

Number of possibilities: 512

PMT1.PRO 12 Ser Tyr Tyr Gln Tyr Glu Gly Ala 19

Oligonucleotide sequence: ntaytaycartaygargngc

2023141744

PMT1. PRO

VED-8

-- PAGE 2

Oligonucleotide sequence: taytaycartaygarggngcn

PROBE#: 11

Melting Temp: 54

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: aytaycartaygarggngcnt

PROBE#: 12

Melting Temp: 54

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: ytaycartaygarggngcntt

PROBE#: 13

Melting Temp: 54

Number of possibilities: 512

PMT1.PRO 14 Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: taycartaygarggngcntty

PROBE#: 14

Melting Temp: 52

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala 19

Oligonucleotide sequence: aytaycartaygarggngcn

PROBE#: 15

Melting Temp: 52

Number of possibilities: 512

PMT1.PRO 13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: ytaycartaygarggngcnt

PROBE#: 16

Melting Temp: 52

Number of possibilities: 512

PMT1.PRO 14 Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: aycartaygarggngcntty

PROBE#: 17

Melting Temp: 52

Number of possibilities: 512

PMT1.PRO 15 Gln Tyr Glu Gly Ala Phe Leu 21

Oligonucleotide sequence: cartaygarggngcnttyt

PROBE#: 18

Melting Temp: 50

Number of possibilities: 512

PMT1.PRO 12 Ser Tyr Tyr Gln Tyr Glu Gly 18

Oligonucleotide sequence: wantaytaycartaygargg

PROBE#: 19

Melting Temp: 50

Number of possibilities: 512

PMT1.PRO 12 Ser Tyr Tyr Gln Tyr Glu Gly Ala 19

Oligonucleotide sequence: ntaytaycartaygarggng

2023141745

PMT 1. PRO

VED - 9

-- PAGE 3

Oligonucleotide sequence: ytaycartaygarggngcn

PROBE#: 21

Melting Temp: 50

Number of possibilities: 512

PMT1.PRO

14 Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: ycartaygarggngcntty

PROBE#: 22

Melting Temp: 50

Number of possibilities: 512

PMT1.PRO

15 Gln Tyr Glu Gly Ala Phe Leu 21

Oligonucleotide sequence: cartaygarggngcntty

PROBE#: 23

Melting Temp: 60

Number of possibilities: 1024

PMT1.PRO

13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20 *

Oligonucleotide sequence: taytaycartaygarggngcntty

PROBE#: 24

Melting Temp: 58

Number of possibilities: 1024

PMT1.PRO

12 Ser Tyr Tyr Gln Tyr Glu Gly Ala 19

Oligonucleotide sequence: sntaytaycartaygarggngc

PROBE#: 25

Melting Temp: 58

Number of possibilities: 1024

PMT1.PRO

13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: aytaycartaygarggngcntty

PROBE#: 26

Melting Temp: 58

Number of possibilities: 1024

PMT1.PRO

14 Tyr Gln Tyr Glu Gly Ala Phe Leu 21

Oligonucleotide sequence: taycartaygarggngcnttyt

PROBE#: 27

Melting Temp: 56

Number of possibilities: 1024

PMT1.PRO

13 Tyr Tyr Gln Tyr Glu Gly Ala Phe 20

Oligonucleotide sequence: ytaycartaygarggngcntty

PROBE#: 28

Melting Temp: 56

Number of possibilities: 1024

PMT1.PRO

14 Tyr Gln Tyr Glu Gly Ala Phe Leu 21

Oligonucleotide sequence: taycartaygarggngcnttyt

PROBE#: 29

Melting Temp: 56

Number of possibilities: 1024

PMT1.PRO

14 Tyr Gln Tyr Glu Gly Ala Phe Leu 21

Oligonucleotide sequence: aycartaygarggngcnttyt

2023141746

PMT2.PRO



-- PAGE 1
PROBE#: 1
Melting Temp: 40
Number of possibilities: 768
PMT2.PRO 6 Val Ala Ile Ser Gly 10 ✓
Oligonucleotide sequence: gtngcnathwsngg

PROBE#: 2
Melting Temp: 38
Number of possibilities: 768
PMT2.PRO 7 Ala Ile Ser Gly Tyr 11
Oligonucleotide sequence: gcnathwsnggnta

PROBE#: 3
Melting Temp: 38
Number of possibilities: 1024
PMT2.PRO 9 Ser Gly Tyr Arg Asp 13
Oligonucleotide sequence: wsnggntaymgnga

PROBE#: 4
Melting Temp: 38
Number of possibilities: 1024
PMT2.PRO 9 Ser Gly Tyr Arg Asp 13
Oligonucleotide sequence: snggntaymgngay

2023141747

PMT3.PRO

VED - II

-- PAGE 1

PROBE#: 1

Melting Temp: 46

Number of possibilities: 192

PMT3.PRO 1 Ala Asp Ile Glu His Tyr 6

Oligonucleotide sequence: gcngayathgarcaytay

PROBE#: 2

Melting Temp: 48

Number of possibilities: 384

PMT3.PRO 1 Ala Asp Ile Glu His Tyr Ser 7

Oligonucleotide sequence: gcngayathgarcaytayw

PROBE#: 3

Melting Temp: 46

Number of possibilities: 384

PMT3.PRO 1 Ala Asp Ile Glu His Tyr 6

Oligonucleotide sequence: ngcngayathgarcayta

PROBE#: 4

Melting Temp: 52

Number of possibilities: 768

PMT3.PRO 1 Ala Asp Ile Glu His Tyr Ser 7

Oligonucleotide sequence: gcngayathgarcaytayws

PROBE#: 5

Melting Temp: 48

Number of possibilities: 768

PMT3.PRO 1 Ala Asp Ile Glu His Tyr 6

Oligonucleotide sequence: ngcngayathgarcaytay

PROBE#: 6

Melting Temp: 48

Number of possibilities: 768

PMT3.PRO 1 Ala Asp Ile Glu His Tyr Ser 7

Oligonucleotide sequence: cngayathgarcaytayws

PROBE#: 7

Melting Temp: 48

Number of possibilities: 768

PMT3.PRO 2 Asp Ile Glu His Tyr Ser Lys 8

Oligonucleotide sequence: gayathgarcaytaywsnaa

PROBE#: 8

Melting Temp: 46

Number of possibilities: 768

PMT3.PRO 2 Asp Ile Glu His Tyr Ser Lys 8

Oligonucleotide sequence: gayathgarcaytaywsna

2023141748

PMT 4. PRO

VED - 12

-- PAGE 1

PROBE#: 1

Melting Temp: 52

Number of possibilities: 384 ✓

PMT4.PRO 2 Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: athacngaraayggnttygc

PROBE#: 2

Melting Temp: 50

Number of possibilities: 384

PMT4.PRO 2 Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: thacngaraayggnttygc

✓ PROBE#: 3

Melting Temp: 56

Number of possibilities: 512

PMT4.PRO 3 Thr Glu Asn Gly Phe Ala Gly 9

Oligonucleotide sequence: acngaraayggnttygcngg

PROBE#: 4

Melting Temp: 54

Number of possibilities: 512

PMT4.PRO 3 Thr Glu Asn Gly Phe Ala Gly 9

Oligonucleotide sequence: cngaraayggnttygcngg

PROBE#: 5

Melting Temp: 52

Number of possibilities: 512

PMT4.PRO 3 Thr Glu Asn Gly Phe Ala Gly 9

Oligonucleotide sequence: acngaraayggnttygcng

PROBE#: 6

Melting Temp: 50

Number of possibilities: 512

PMT4.PRO 3 Thr Glu Asn Gly Phe Ala Gly 9

Oligonucleotide sequence: cngaraayggnttygcng

PROBE#: 7

Melting Temp: 50

Number of possibilities: 512

PMT4.PRO 3 Thr Glu Asn Gly Phe Ala Gly 9

Oligonucleotide sequence: ngaraayggnttygcngg

PROBE#: 8

Melting Temp: 50

Number of possibilities: 512

PMT4.PRO 4 Glu Asn Gly Phe Ala Gly 9

Oligonucleotide sequence: garaayggnttygcngg

✓ PROBE#: 9

Melting Temp: 58

Number of possibilities: 768

PMT4.PRO 1 Phe Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: ttyathacngaraayggnttygc

PROBE#: 10

Melting Temp: 56

Number of possibilities: 768

PMT4.PRO 1 Phe Ile Thr Glu Asn Gly Phe Ala 8

2023141749

PMT 4. PRO

VED-13

-- PAGE 2

Oligonucleotide sequence: tyathacngaraayggnttygc

PROBE#: 11

Melting Temp: 54

Number of possibilities: 768

PMT4.PRO

1 Phe Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: ttyathacngaraayggnttyg

PROBE#: 12

Melting Temp: 54

Number of possibilities: 768

PMT4.PRO

1 Phe Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: yathacngaraayggnttygc

PROBE#: 13

Melting Temp: 52

Number of possibilities: 768

PMT4.PRO

1 Phe Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: tyathacngaraayggnttyg

PROBE#: 14

Melting Temp: 50

Number of possibilities: 768

PMT4.PRO

1 Phe Ile Thr Glu Asn Gly Phe 7

Oligonucleotide sequence: ttyathacngaraayggntty

PROBE#: 15

Melting Temp: 50

Number of possibilities: 768

PMT4.PRO

1 Phe Ile Thr Glu Asn Gly Phe Ala 8

Oligonucleotide sequence: yathacngaraayggnttyg

PROBE#: 16

Melting Temp: 56

Number of possibilities: 1024

PMT4.PRO

4 Glu Asn Gly Phe Ala Gly Arg 10

Oligonucleotide sequence: garaayggnttygcnggnmg

PROBE#: 17

Melting Temp: 52

Number of possibilities: 1024

PMT4.PRO

4 Glu Asn Gly Phe Ala Gly Arg 10

Oligonucleotide sequence: garaayggnttygcnggnm

PROBE#: 18

Melting Temp: 52

Number of possibilities: 1024

PMT4.PRO

4 Glu Asn Gly Phe Ala Gly Arg 10

Oligonucleotide sequence: araayggnttygcnggnmg

PROBE#: 19

Melting Temp: 50

Number of possibilities: 1024

PMT4.PRO

4 Glu Asn Gly Phe Ala Gly Arg 10

Oligonucleotide sequence: raayggnttygcnggnmg

2023141750

REVERSE PRIMERS=c mRNA EcoR1?

PMT4.PRO (15 kd CNBr fragment)

1	2	3	4	5	6	7	8	9
Met	Phe	Ile	Thr	glu	Asn	Gly	Phe	Ala
10	11	12	13	14				
GLY	Arg	Ser	Gly	Arg				

1, 5'-C CIG CRA AIC CRT TYT CIG TRA
TRA ACA T -3',10(P-15).

1, 5'-ATC GAATTC CIG CRA AIC CRT TYT
CIG TRA TRA ACA T -3,'10(P-16).

5'-ATC GAATT CCR TTY TCN GTD ATR
AA CAT-3'(P-21, REVERSE OF P-17)

5'-AAN CCR TTY TCN GTD ATR AAC AT-
3'(P-22, REVERSE OF P-18).

2023141751

FORWARD PRIMERS=mRNA BamH1?
 PMT4.PRO (15 kd CNBr fragment)

R:AG,Y:CT,H:ATC,D:ATG,B:GTC,N:AGCT

1	2	3	4	5	6	7	8	9
Met	Phe	Ile	Thr	glu	Asn	Gly	Phe	Ala
10	11	12	13	14				
GLY	Arg	Ser	Gly	Arg				

1 AAC GGATCC ATG TTY ATH ACN GAR
 AAY GG 7(P-17)

MELTING TEMPERATURE: 50

Number of possibilities: 96

1 ATG TTY ATH ACN GAR AAY GGN TT
 8(P-18)

Melting Temperature: 56

Number of Possibilities: 384

FORWARD PRIMERS=mRNA,BamH1?

PMT1.PRO(N-29 mer OF 60 kd PROTEIN)

LSSNFLFGTASSYYQYEG AFL
SDGVGLSN

13 Tyr Tyr Gln Tyr Glu Gly Ala 19
ATC GGA TCC TAY TAY CAR TAY GAR
GGI GC (P-19)

Melting temp: 52

Number of possibilities: 32

13 TAY TAY CAR TAY GAR GGI GCI TT
20.(P-20)


Melting temperature: 58

Number of possibilities: 32


ONE-SIDED PCR

Cloning of the 3' End of a cDNA with 5' Sequence Information

I. REVERSE TRANSCRIPTASE REACTION




II. AMPLIFICATION BY PCR



2023141754

Anchor PCR





Loh, E.
1990

AMERICAN UNIVERSITY

2023141755

M. Frohman, PNAS (1988) 85:8998-9002



2023141756

-40 -30 -20 -10 -1 1
 CGTGTGTCCTCAAAGTGGCCACAGGGCGCTGACAGCAGC
 10 20 30 40 50 60 70 80 90 100 110 120 130 140 150 160 170 180 190 200 210 220 230 240 250 260 270 280 284
 CAG GCG GCG GCA GTG CCC GAC TCA GAC CCG GGC CTG GCG GCG GGC ACA GAC CCG AGT
 Gln Ala Ala Gly Ala Val Pro Asp Ser Asp Pro Gly Leu Ala Ala Val Ser Ser Ala Tyr
 CAG CGC TTT GAG CCC GCG CTC TAC CTC CGC AAC AAC TAC GCG CCC CCG AGG GGG GAC CTG
 Gln Arg Phe Glu Pro Arg Ala Tyr Leu Arg Asn Asn Tyr Ala Pro Pro Arg Gly Asp Leu
 AGC TGC CCC GAC GGC GTC GGG CCT TGG AAG CTG CCG TGC TTT GCT CAG ACC TTC GCC ACC
 Ser Cys Pro Asp Gly Val Gly Pro Trp Lys Leu Arg Cys Leu Ala Gln Thr Phe Ala Thr
 GGT GAG GTG TCT GGC CGC ACC CTC ATT GAC ATC GGT TCA GGA CCC ACT ATA TAC CAG CTG
 Gly Glu Val Ser Gly Arg Thr Leu Ile Asp Ile Gly Ser Ser Gly Pro Thr Ile Tyr Gln Leu
 CTC AGC GCC TGT GCC CAC TTT GAG GAC ATC ACC ATG ACA GAT TTC CTG GAG GTG AAC CGC
 Leu Ser Ala Cys Ala His Phe Glu Asp Ile Thr MET Thr Asp Phe Leu Glu Val Asn Arg
 CAG GAG CTG AGG CTC TGG CTG CGA GAA GAG CCT GGG GCT TTC GAC TGG AGC GTG TAC AGC
 Gln Glu Leu Arg Leu Trp Leu Arg Glu Glu Pro Gly Ala Phe Asp Trp Ser Val Tyr Ser
 CAG CAT GTC TGC CTC ATC GAG GGC AAG GGG GAA TCC TGG CAG GAG AAG GAG TGC CAG CTG
 Gln His Val Cys Leu Ile Glu Gly Lys Gly Glu Ser Trp Trp Gln Glu Lys Glu Cys Gln Leu
 CGA GCC AGG GTG AAG AGG ATC CTG CCC ATC GAT GTG CAC CCG CCC CAG CCC CTG GGT GCT
 Arg Ala Arg Val Lys Arg Ile Leu Pro Ile Asp Val His Arg Pro Gln Pro Leu Gly Ala
 GGA GGC CTG GCA CCC CTG CCT GCC GAC GCC CTG GTC TCT GCC TTC TGC CTG GAG GCT GTG
 Gly Gly Leu Ala Pro Leu Pro Ala Asp Ala Leu Val Ser Ala Phe Cys Leu Glu Ala Val
 AGT CCA GAC CTG GCC AGC TTC CAG CGG GCC CTG GAC CAC ATC ACC ACA CTG CTG AGG CCT
 Ser Pro Asp Leu Ala Ser Phe Gln Arg Ala Leu Asp His Ile Thr Thr Leu Leu Arg Pro
 GGG GGG CAC CTC CTC ATC GGA GCC CTG GAG GAG TCA TGG TAC CTG GCT GGG GAG GCC
 Gly Gly His Leu Leu Ile Gly Ala Leu Glu Glu Ser Trp Tyr Leu Ala Gly Glu Ala
 200 210 220 230 240 250 260 270 280 284
 AGG CTG GCG GTG CCC GTG CCG GAG GAG GTG AGG GAG GCG GCG CTG GTG CCG AGC GCT
 Arg Leu Ala Val Val Pro Val Arg Glu Glu Val Arg Glu Ala Leu Val Arg Ser Ala
 ATG AGG TGC GGG ATC TGC GCA CTA CAC CAT CAC GCG GCT TGC CCA CCT TCA GAC AGG TGT
 MET Arg Cys Gly Ile Cys Ala Leu His Leu His His Ala Cys Pro Pro Ser Asp Arg Cys
 280 284
 AGA CGA TGT CAA GGG CAT CTT CAC CTC GCC CAG AAG GAG GTG GGG GTG TGA GCCCC
 Arg Arg Cys Gln Gly His Leu Leu His Leu Ala Gln Lys Val Gly Val Ter




FIG. 6. Comparison of the location of regions of sequence similarity of the human D-Asp/L-isoAsp and other methyltransferases along their polypeptide chains. The homology regions, connected by vertical lines, are: Region I (diagonal shading), Region II (dotted shading), and Region III (diamond shading). Sequences are detailed in Tables IX (Region I) and X (Regions II and III). Although not listed, Region I is located at a similar position in the other bacterial RNA adenine methyltransferases.

ACTION PLAN 1992

THESE ARE ONLY SUGGESTED PLANS.
PRIORITY DECREASES WITH THE
INCREASING ASSIGNED NUMBER

1. Obtain aminoacid sequence (~~c-terminal~~)
of PMT.

2. Use aminoacid sequence information to
synthesize primers for use in PCR.

3. USE PCR to isolate desired sequences
from the tobacco root RNA and DNA.

a. RNA PCR(Perkin-Elmer kit, Forward
and Reverse primers)

b. DNA pcr(Perkin-Elmer kit, Forward
and Reverse Primers)

c. ANCHORED PCR:AMPLIFICATION
WITH SINGLE-SIDED SPECIFICITY

1. Take tobacco root poly A
RNA. Synthesise complementary strand of

2023141759

DNA using oligo dT as primer .Destroy RNA by alkali hydrolysis.Amplify Using combinatins of forward primers and oligo dT25 .CLONE the PCR PRODUCT.SEQUence the insert. Enter the sequence into computer.use DNA STAR(findpro programme) to find PMT1.PRO,PMT2.PRO, PMT3.PRO PMT4. PRO sequences.

2.ANCHORED PCR of cDNA with OLIGO(dG)tails.Synthesise cDNA using REVERSE PCR PRIMERS or oligo dT.Add oligo dG to the newly synthesised c DNA strand.Amplify using REVERSE PCR PRIMERS and OLIGO dC.Clone the PCR product AND analyse as in 1.

d.LIGATION-MEDIATED ANCHORED PCR.

Digest the genomic DNA with a Mbo1 OR its isoschizomer Sau3A.

If possible, Fill one base to prevent self-Ligation.

Ligate an anchor-primer on the ends of the DNA fragments.

Amplify with a specific primer as well as the anchor.

The choice of which strand the primer is derived from determines which side of the fragment is amplified .

This can be done for several enzymes that leave cohesive ends with a four-base 5'overhang. This could allow isolation of the promoter for the gene of interest.

4. GENERATION OF SINGLE-STRANDED LABELED PRIMERS AND THEIR USE IN NORTHERN AND SOUTHERN HYBRIDIZATION

5. IMMUNOSCREENING OF EXPRESSION LIBRARY >

2023141761

SEQUENCE INFORMATION

1. N-terminus of intact 60 kD protein, electroblot/sequenced. (code:PNMT1.PRO)

367

PMT L S S N F L F G T A S S Y Q Y E G A F L S D G V G L S M --- COOH
 Deg 6 6 6 2 2 6 2 4 4 4 6 6 2 2 2 2 4 4 2 6 6 2 4 4 4 6 6 2

Alignments:

PNMT (PET) S(16) S S Y Q
 PNMT (Phosphatidyl) F Y G D F F L S

2. 6.2 kD CNBr Fragment of 60 kD protein, electroblot/sequenced (F-SG-T MOTIF, region I), (a met would precede all of the CNBr fragments). (code:PNMT2.PRO)

PMT M E P X E V A I S G Y R D X I --- COOH
 Deg 2 2 4 2 4 4 3 6 4 2 6 2 4

Alignments:

PNMT C/G I E E I G S G K G H F I
 PNMT (65) E A T G E S G R I
 PNMT (75) L I D I G S G P T Y N L

3. 3 kD CNBr Fragment of 60 kD after removal of N-terminus (DAL - Region II). (code:PNMT3.PRO)

PMT A D I E H S K L I D A L X I K G I Q F --- COOH
 Deg G L L A P 2 6 2 6 3 2 4 6 3 2 4 3

Alignments:

PNMT (125) X S K G
 PNMT (174) P A D A L
 AdoHcyHMT K D A I

4. 15 kD CNBr Fragment of 60 kD protein. (code:PNMT4.PRO)

PMT F I T E N G F A G R S G R P (?) --- COOH
 Deg 2 3 4 2 2 4 2 4 4 6 6 4 6

5. Possible sequence order of fragments obtained for 60 kD protein.

From this speculative "motif" examination, a deduced composite can be drawn for putative PMT as follows (cf. PNMT, PNAS 83,5455 (1986):

N-terminus (~3 kD) ----- 6.2 kD fragment ----- other 3 kD fragment ----- 15 kD fragment ---- COOH

2023141762

IN ROOTING STAGE

172 TOTAL

121	VST
7A	46A
1A	49A
13A	19A
47A	9A
50A	17.3A
10.6A	3A & 10A

2023141763

IN ROOTING STAGE

172 TOTAL

121	VST
7A	46A
1A	49A
13A	19A
47A	9A
50A	17.3A
10.6A	3A & 10A

2023141764

GENERATION OF TRANSGENIC PLANTS
CURRENT STATUS

IN GREENHOUSE

77 TOTAL

121	13A
VST	46A
17.2A	1A
49A	12A
*47A	*19A
*3A & 10A	

2023141765

After 10-14 days culture, separate developing calli and transfer to fresh medium



As soon as morphologically acceptable shoots appear, transfer shoots to rooting medium containing 500 mg/L carbenicillin and 100 mg/L kanamycin



After rooting, transgenic plantlets can be transferred to greenhouse or maintained as meristem cultures indefinitely until greenhouse space is available

LEAF DISC TRANSFORMATION

Grow transformed bacterial clones
overnight in suspension culture



Excise leaf material from axenic Burley 21
plantlets and cut lamina into 5-10 mm
sections



Dip sections into cell culture, blot on sterile
gauze and transfer to shoot induction
medium without antibiotics



After 2 days incubation (28⁰ C, 16 hr
photoperiod, 5000 lux), rinse sections in
500 mg/L carbenicillin, blot and transfer to
shooting medium with 500 mg/L
carbenicillin and 300 mg/L kanamycin

2023141767

Electroporate cells at 100 μ F capacitance,
250 μ sec discharge interval, 480 v, .3mm
gap for final field strength of 16000 v/cm



Dilute cells 10-fold with nonselective culture
medium and incubate at rt for 30-60
minutes



Plate 10 μ l and 70 μ l on Luria + 100 μ g/ml
kanamycin



After incubating for 2-3 days at 28⁰ C,
select 3 isolates from each construct and
generate pure cultures for leaf
transformation

AGROBACTERIUM ELECTROPORATION

Grow overnight suspension culture of
Agrobacterium tumefaciens str. LBA4404



Wash cells 2x in 1mM HEPES + 10%
glycerol with a final wash in 10% glycerol
alone (reduce volume to 500 μ l)



Mix 2 μ l of DNA with 25 μ l cells



Transfer 8 μ l of mix to Hoefer bacterial
electrode (5.6mm)

GENERATION OF TRANSGENIC PLANTS
VIA
PLANT TISSUE CULTURE

2023141770

GENERATION OF TRANSGENIC PLANTS
VIA
PLANT TISSUE CULTURE

ASK AND YE
SHALL RECEIVE

BBM
→

2023141771

ASK AND YE
SHALL RECEIVE

BGM
→

2023141772

**GENERATION OF TRANSGENIC PLANTS
VIA
PLANT TISSUE CULTURE**

2023141773